



# VUETRACK RF

## What Is VueTrack-RF?


VueTrack-RF™ is a groundbreaking solution that tracks medical supplies, assets and people in real time with 99.5% accuracy, using advanced ultra-high frequency (UHF) RFID antennas and inexpensive, versatile, passive tags. VueTrack-RF turns any space – whether a supply room, procedure room, or warehouse – into a fully controlled inventory space using existing shelving. VueTrack-RF merges VueTrack's™ proven technology for tracking, documenting and reporting inventory transactions and data with UHF RFID, another leading and well-recognized data capture and tracing technology. Together, these two technologies are bringing transformations to healthcare that other industries have already embraced.

## How Does VueTrack-RF Work?

1. Employs zonal and steerable “smart” (BESPA) antennas along with handheld RFID readers with VueTrack-RF location processors to continuously and automatically track the movement and location of all clinical supplies.
2. Uses safe and inexpensive global standards compliant (GS1 EPC) passive Gen 2 RFID tags from any manufacturer. These tags meet UDI compliance regulations when used with our RFID Registration and Encoding Solution.
3. Works with the Kanban system: RFID-tagged Kanban cards are placed at PAR locations during stocking and then dropped into a designated drop box to trigger the reordering of that item. This is particularly effective for tracking small or inexpensive med surg or commodity items.

## Why Use VueTrack-RF?

- Makes the monitoring and accounting of clinical supplies efficient, accurate and easy. Hands-free, no human interaction, virtually error-proof.
- Alerts users to a product's location as well as unapproved or unexpected movements. Tracks assets, patients and personnel with equal agility.
- No proprietary RFID storage cabinets needed. No costs or hassles associated with bulky cabinets, existing storage to dispose of, or departments to remodel. Small antennas are strategically placed on walls or in the ceilings; everything else remains the same.
- Supports UDI regulatory compliance and meets UDI requirements for all classes of devices when our RFID Registration and Encoding Solution is implemented with VueTrack-RF.
- Ensures product availability, optimizes par levels, and reduces inventory size.
- Increases billed revenue.
- Increases data accuracy, thereby immediately optimizing inventory and replenishment.
- Automates expiration tracking to significantly decrease the occurrence of products expiring on shelves.
- Ensures a quick response to recalls.
- ROI realization begins with the first procedure. VueTrack-RF pays for itself several times over within the first year. Benefits continue to accrue thereafter.



## Benefits for Manufacturers and the Supply Chain

VueTrack-RF automatically tracks Work-In-Progress status, inventory and movement throughout manufacturing facilities. VueTrack-RF enables manufacturers to:

- ✓ Track their own manufacturing process/supply chain at all stages, like leaders in other industries.
- ✓ Trace products hands-free through the warehouse to their final destination.
- ✓ Improve forecasting and create reliable prediction models for production and distribution.
- ✓ Gain real-time knowledge of the whereabouts of their products and what quantities are at which hospitals when using VUEMED's Vendor Managed Inventory (VMI) program, thereby enabling manufacturers to:
  - Control expirations and eliminate waste.
  - Set appropriate SKU-based par levels to ensure product availability.
  - Customize composition of consigned and non-consigned inventory.
  - Manage any recalls swiftly and efficiently.
  - Track products all the way to the patient.
- ✓ Get exceptional usage data which can improve offerings and services to customers and generate substantial savings and market share growth.



Connect with us at [www.vuemed.com](http://www.vuemed.com)

