

VUEBIN

What is VueBin?

VueBin™ is a Cloud-based RAIN RFID-powered Kanban inventory management system. It allows users to manage all Med-Surg or typically high velocity and inexpensive or non-chargeable supplies that are essential to all procedures or case carts. VueBin elevates the Kanban inventory control system to its highest ever degree of automation thanks to real-time information which generates replenishment requests and item reorders with unparalleled accuracy and diligence.

How does VueBin work?

At the point of storage, VueBin uses VUEMED's unique and proprietary RFID switches (the "VueSwitch") to indicate whether one or both bins in a 2-bin supply management approach needs replenishment. All the user has to do is slide the RFID switch with one finger when a bin becomes empty and immediately an RFID signal for replenishment is automatically sent to VueBin which alerts staff of the exact location of the SKU to replenish, down to its specific rack and exact bin on the shelf, as well as the par quantity that the bin needs to be replenished to. Because of the automation, accuracy and immediacy of the replenishment signal, you now have peace of mind knowing that stock-outs are virtually eliminated.

As part of VUEMED's ecosystem of RAIN RFID technologies, VueBin can work either as a stand-alone Kanban system for your organization or together with our other technologies for managing supplies, implants and medical devices, surgical instruments, and orthopedic, trauma and spine implants and trays.


Why use VueBin?

VueBin reports and aggregates instantly all replenishment needs with the most accurate quantities. VueBin then automatically decreases such quantities of SKUs from the general stock or warehouse supplies and generates an order to purchase more as needed.

VueBin is built with a visualizer that helps see by department, by room, and by rack the status of each bin. It also collects and analyzes data about the replenishment velocity, par quantities, and usage profile of each location to help make better decisions about how to optimize the Kanban system and anticipate the needs of users.

Let VueBin help you achieve:

1. Automated bin replenishment system
2. Accurate inventory control
3. Full visibility of all bins and in real time
4. An end to stock-outs
5. Powerful data analytics and reports
6. Lean and optimized supplies management



Shelf	Product	Max Qty	Replenish # of Bins	Qty	Velocity
Shelf 1	Hypodermic Syringe w/out Needle, 30 cc, Luer Slip Tip Sufe 55-395	40	0/2	0	0d 16h 6m
Shelf 2	Caesal Springs, 10.2 cm, 10.2 cm Sufe 16-0242	30	1/2	15	0d 16h 21m
Shelf 4	Hypodermic Medical Taps, 1 in, 10 yd, paper Sufe 15-30-1	20	2/2	20	0d 16h 21m
Shelf 5	CLINEL 2 in, 75 in, bulk Sufe 2242	40	1/2	20	0d 0h 20m
Shelf 7	ALD, 5.0 in, 0.0 in, Sterile Sufe N0N214504	40	2/2	40	0d 16h 21m
Shelf 8	CLINEL 4 in, 4 in, 16 ply Sufe 2733	30	2/2	30	0d 16h 20m